

HistoQu (plugin for QuPath) installation guide





- 1) Download the HistoJSDK .jar from <https://host.pathomation.com/plugins/qupath/>

[host.pathomation.com - /plugins/qupath/](https://host.pathomation.com/plugins/qupath/)


[\[To Parent Directory\]](#)

11/20/2020 1:58 PM	11204967	HistoJSDK .jar	←
11/23/2020 1:02 PM	488922394	QuPath-0.2.3.zip	
11/23/2020 12:42 PM	168	web.config	


- 2) Open the Qupath Application file

	QuPath-0.2.3 (console).exe	11/09/2020 9:58	Application	115 KB
	QuPath-0.2.3 (console).ico	11/09/2020 9:58	Icon	98 KB
	QuPath-0.2.3.exe	11/09/2020 9:57	Application	114 KB
	QuPath-0.2.3.ico	11/09/2020 9:57	Icon	98 KB

- 3) Drag and drop the plugin into the application window of QuPath



- 4) QuPath automatically adds the Plugin to its extension menu



5) **Go to <https://adoptopenjdk.net/archive.html?variant=openjdk14&jvmVariant=hotspot>**



6) **Download JRE 14.0.2 +12 ZIP file**

macOS x64	JDK	JDK (195 MB)	Checksum
	JRE	JRE (49 MB)	Checksum
Windows x86	JDK	JDK (184 MB)	Checksum
	JRE	JRE (48 MB)	Checksum
Windows x64	JDK	JDK (206 MB)	Checksum
	JRE	JRE (53 MB)	Checksum


7) **Open the Runtime folder from Qupath**

Name	Date modified	Type	Size
app	11/09/2020 9:58	File folder	
runtime	18/02/2021 10:06	File folder	
.jpackage.xml		XML Document	1 KB
api-ms-win-core-console-l1-1-0.dll	11/09/2020 9:57	Application extens...	19 KB
api-ms-win-core-datetime-l1-1-0.dll	11/09/2020 9:57	Application extens...	19 KB
api-ms-win-core-debug-l1-1-0.dll	11/09/2020 9:57	Application extens...	19 KB
api-ms-win-core-errorhandling-l1-1-0.dll	11/09/2020 9:57	Application extens...	19 KB
api-ms-win-core-file-l1-1-0.dll	11/09/2020 9:57	Application extens...	22 KB
api-ms-win-core-file-l1-2-0.dll	11/09/2020 9:57	Application extens...	19 KB
api-ms-win-core-file-l2-1-0.dll	11/09/2020 9:57	Application extens...	19 KB


Tooltip for api-ms-win-core-date-time-l1-1-0.dll:
Date created: 23/11/2020 13:01
Size: 321 MB
Folders: bin, conf, include, jmods, legal, lib
Files: release

8) **Delete all files from this folder**

9) **Extract the JRE 14.0.2 +12 ZIP file to a temporary folder**



10) **Copy the content of the "jdk-14.0.2+12-jre" folder to your empty Runtime folder in Qupath**



11) **Now, you're good to go!**

